

# SHARPS BOXES

## Assembly and disposal

**Lock**, by twisting the top until the red lug clicks into position.

The lid **MUST** be clicked into place on all sides during assembly.

The identifying label **MUST** be completed at the start of use/assembly date and disposal date.



**These are Legal Requirements!**



**Lock**, by pressing the flap down to click, followed by pressing the D-shape tab into the locked position.

Dispose of sharps containers once in use for three months or when 'fill line' is reached.

Put station ID sticker here.

The **station ID** and date sticker must be attached **IN ADDITION** to the completed box details (circled above). **NOT INSTEAD OF!**